
## **COPERNIC RC** *Venous Remodeling Balloon*





## **COPERNIC RC** Venous Remodeling Balloon

- Indication: The "Venous remodeling technique" in the treatment of dural arteriovenous fistulas was created by Pr. René Chapot (Alfried Krupp Krankenhaus, Essen, Germany) to preserve the sinus during embolization of the meningeal feeders.
- At first, this technique was done by placing angioplasty PTA balloons into the sinus, which were very difficult to position and not compliant enough. To overcome this problem, BALT was asked to create a special venous balloon: the COPERNIC RC. With a long length and a large diameter, this balloon is perfect for this procedure thanks to its **ease of navigation** and its **high suppleness and compliance**.
- The COPERNIC RC is introduced via a femoral venous puncture into the sinus by venous navigation with a 6F guiding catheter FARGOMAX. The injection of the embolization product is made by arterial approach. Inflation of the balloon when the embolization product reaches the vein induces a perfect rerouting of the product to the other feeders.



- Occlusion of sinus with the balloon's long length
- Easy navigation thanks to HYDROSPEED coating on the balloon
- Easy control of the procedure thanks to its high visibility
- DMS0 compatibility



	Reference	Catheter length (cm)	Balloon length (mm)	Balloon diameter (mm)	maximum outer diameter (mm)	Guidwire	Compatible guiding catheter
N.	COPERNIC8X80RC	160	80	8	4,35F (1.45mm)	Compatible with Transend14 guidwire Stryker*	FARGO MAX 6F (I.D. : 1.78mm / .070" )
	COPERNIC10X80RC	160	80	10	5F (1,6mm)		

\* Other guidewires may induce risk of leakage.



10, rue de la Croix Vigneron 95160 Montmorency • France Tél. : +33 (0)1 39 89 46 41 Fax : +33 (0)1 34 17 03 46 www.balt.fr