

# ordering information

Developed portfolio with unique references

## Valvuloplasty

Reference	Balloon Ø x L	Lc*	Ø IVA	GW**
CBV12x30/110	12 x 30	110	7F	.038"
CBV15x30/110	15 x 30	110	8F	.038"
CBV15x40/110	15 x 40	110	8F	.038"
CBV16x40/110	16 x 40	110	8F	.038"
CBV18x30/110	18 x 30	110	8F	.038"
CBV18x40/110	18 x 40	110	8F	.038"
CBV18x45/110	18 x 45	110	8F	.038"
CBV20x30/110	20 x 30	110	8F	.038"
CBV20x40/110	20 x 40	110	8F	.038"
CBV20x45/110	20 x 45	110	8F	.038"
CBV23x45/110	23 x 45	110	10F	.038"
CBV25x40/110	25 x 40	110	10F	.038"
CBV25x50/110	25 x 50	110	10F	.038"
CBV28x50/110	28 x 50	110	12F	.038"
CBV30x60/110	30 x 60	110	12F	.038"
CBV35x60/110	35 x 60	110	15F	.038"
CBV40x40/110	40 x 40	110	15F	.038"

\* Lc: Length of the catheter

\*\* GW: Guidewire

## Valvuloplasty small sizes

Reference	Balloon Ø x L	Lc*	Ø IVA	GW**
CBVP8x20/70	8 x 20	70	6F	.035"
CBVP10x20/70	10 x 20	70	6F	.035"
CBVP12x20/70	12 x 20	70	7F	.035"
CBVP15x25/70	15 x 25	70	8F	.035"

## Peripheral Angioplasty

Reference	Balloon Ø x L	Lc*	Ø IVA	GW**
CBP8x30/110	8 x 30	110	6F	.035"
CBP10x30/110	10 x 30	110	6F	.035"
CBP12x40/110	12 x 40	110	7F	.035"


Cristal balloons are catheters designed for peripheral angioplasty and valvuloplasty. Class III CE0297 in compliance with Medical Device Directive (MDD 93/42/EEC amended by 2007/47/EC). Manufactured by BALT Extrusion SAS. Carefully read the instructions for use before use.  
First CE marking: 1998. The content of this document, in particular data, information, trademarks and logos is BALT SASs and affiliates' sole property. © 2019 BALT SAS and affiliates, all rights reserved. All representation and/or reproduction, whether in part or in full, is forbidden and would be considered a violation of BALT SAS and affiliates' copyrights and other intellectual proprietary rights.  
This document with associated pictures is non-contractual and is solely dedicated to healthcare professionals and BALT's distributors (BALT's supplier's distributors). The products commercialized by BALT shall exclusively be used in accordance with the instructions for use included in the boxes. DC006GB (03/2020)

Balt  
10, rue de la Croix Vigneron, 95160 Montmorency France  
Tél. : +33 (0)1 39 89 46 41  
Fax : +33 (0)1 34 17 03 46  
[www.balt.fr](http://www.balt.fr)




# cristal balloon

angioplasty and  
valvuloplasty balloon



# cristal balloon

Designed for peripheral angioplasty & valvuloplasty



High RBP\*

CBV 20mm = 8atm

Over

160,000


cristal balloons used worldwide

A large scale of diameter  
Up to 40mm Ø and 60mm length


Semi-compliant behavior

Controlled conformability  
Around nominal pressure

High RBP\*  
For better reliability against rupture



20x30mm Cristal Balloon Compliance chart



\*Rate burst pressure

peripheral treatment