




CARDIAC INTERVENTION

Worley™ Vein Selector

Advancing
Vessel
Selection




Part of the
Worley Advanced
Telescopic System

Worley™ Vein Selector

PRODUCT FEATURES

- Unique accordion hub design to minimize kinking and enhance guiding capacity
- Proprietary braided shaft design to maximize flexibility
- Soft tapered tip to reduce vessel trauma
- Highly visible under fluoroscopy
- Improve LV lead location.¹



Greater visibility for reliable access to the target vein




Ergonomic wing hub enhances handling and control

ORDERING INFORMATION

Part Number	French Size	Internal Diameter	Length	Curve
57538CS-WOR	5F	.046"/ 1.17 mm	75 cm	Standard
57538CSV-WOR	5F	.046"/ 1.17 mm	75 cm	Vertebral
57538CSHK-WOR	5F	.046"/ 1.17 mm	75 cm	Hook

¹ Jackson KP , Hegland DD, Frazier-Mills C, Piccini JP, Koontz JI, Atwater BD, Daubert JP, Worley SJ. Impact of using a telescoping-support catheter system for left ventricular lead placement on implant success and procedure time of cardiac resynchronization therapy. Pacing Clin Electrophysiol. 2013 May; 36(5):553-8.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.208.4300
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

Merit.com

Austria
0800 295 374

Belgium
0800 72 906 (Dutch)
0800 73 172 (Français)

Denmark
80 88 00 24

Finland
0800 770 586

France
0800 91 60 30

Germany
0800 182 0871

Ireland (Republic)
1800 553 163

Italy
800 897 005

Luxembourg
8002 25 22

Netherlands
0800 022 81 84

Norway
800 11629

Portugal
308 801 034

Russia
+7 495 221 89 02

Spain
+34 911238406

Sweden
020 792 445

Switzerland
(Deutsch)
+41 225180252
(Français)
+41 225948000
(Italiano)
+41 225180035

UK
0800 973 115