Products for Pulmonology



SINGLE USE BOUTIN TROCARS FOR THE PLEURA

For pleural biopsies and punctions

NOVATECH has launched trocars for the pleura in cooperation with Professor C. Boutin (Marseilles, France) in 1998. He redesigned the existing trocar to improve the quality of sampling.

The single use Boutin trocars have a working length of 78 mm. They come sterile and are available in two diameters (2 and 3 mm).

The three-sided mandrel and the improved shape of the hook of the 3 mm trocar give this device a sampling quality higher than with the Abrams needle¹⁾.

2 mm: for pleural punction only3 mm: for pleural biopsy and punction



Boutin Trocar 2 mm: for pleural punction only, with a lateral orifice and two mandrels: sharp resp. blunt.



Boutin Trocar 3 mm: for pleural biopsy and punction, with lateral hook and two mandrels: sharp resp. blunt.



Features

- sterile, ready for use
- convenient dispenser box
- single use, minimizing risk of contamination
- less traumatic due to diameters adapted to intended use
- equipped with a three-way stopcock

REF	Description	Items per box
58023	Disposable BOUTIN Trocar for pleural punction 2 mm, with two mandrels	12
58033	Disposable BOUTIN Trocar for pleural biopsy and punction, 3 mm, with two mandrels	12
STERILE		

www.novatech.fr

page 5

¹⁾ Study by Prof. C. BOUTIN - Hôpital La Conception, Marseilles Presentation: CHEST/Toronto/1998